



Laboratory automation and data management in diagnostics

*qPCR ZIKV detection and
quantification*

PAIN




3 years

Scientific data output
doubles

54%



Of studies cannot be verified




17%

Of scientific data is lost
each year

Research: 1900s vs. 2000s



WE PAPER





Diagnostic setting – early automation

High throughput → first to automate

Automate with software – LIS, LIMS

Automate with instruments



Automation in research?

Generally less throughput → less pressure to automate

Instruments – same as diagnostics

DIGITISATION



LABORATORY DIGITISATION



LABORATORY DIGITISATION

The screenshot displays a Trello board for 'sciNote M&S Meetings'. The board is organized into several columns representing different stages or time periods:

- To-Do:** A list of tasks including Business trip, Outbound marketing, Inbound marketing, Meeting (checked), Support, Sales, Marketing material, and Development.
- Reminders:** A list of reminders such as Newsletter for subscribers, sciNote Blog, Release notes, and a meeting review for 2017.
- Weekly Updates (Weeks 14-17):** Columns showing progress for specific weeks, with cards for tasks like 'Blog articles', 'Webpage improvements', 'Wikipedia', 'Workshop', 'Email campaign', and 'Video infographics'.

The interface includes a top navigation bar with 'Boards', a search icon, the Trello logo, and a user profile for 'Jana Erjavec'. A 'Show Menu' option is visible in the top right corner.

LABORATORY DIGITISATION

Google

Mail

COMPOSE

Inbox (1,574)

Starred

Sent Mail

Drafts (3)

MarketsandMarkets 9:30 am
Environmental Monitoring Market worth \$20.5 Billion by 2020
The emerging markets including China, India, Brazil,
Jure Pecher 9:30 am
Vizualizacija poteka projektov
Živjo Jana, sem šel čez zapiske z najinega pogovora

Environmental Monitoring Market worth \$20.5 Billion by 2020

MarketsandMarkets <Newsletter@marketsandmarkets.com> 9:30 AM (1 minute ago)
to Jana

The emerging markets including China, India, Brazil, Email not displaying correctly?

Add card

1-50 of 9,397



GMail to Trello

offered by nstung@gmail.com

★★★★★ (289) Productivity 29,264 users

Environmental Monitoring Market worth \$20.5 Billion by 2020

MarketsandMarkets <Newsletter@marketsandmarkets.com>

The emerging markets including China, India, Brazil, become attractive destinations for companies engaged in the development and distribution of Environmental Monitoring Products.

Environmental Monitoring Market worth \$20.5 Billion by 2020

Environmental Monitoring Market worth \$20.5 Billion by 2020
Gas, Water, Soil, and Noise Pollution

1-50 of 9,397

janaerjavec | Logout? | Features/Bugs [x] Close

Board.: My Boards » Software development

List: Info To Do Doing Done Temporary on hold

Title: Environmental Monitoring Market worth \$20.5 Billion by 2020

Description: The emerging markets including China, India, Brazil, and Mexico have become attractive destinations for companies engaged in the development and distribution of Environmental Monitoring Products. Email not displaying correctly? View it in your browser. Unsubscribe
Environmental Monitoring Market worth \$20.5 Billion by 2020
The report "Environmental Monitoring Market by Product (Monitors, Sensors, and Software). Application

Link back to Gmail Assign me to this card

Add to Trello card

LABORATORY DIGITISATION



BOOMERANG FOR GMAIL



zapier



Google Sheets



Gmail



Trello



Google Calendar



Slack



MailChimp



Evernote



Twitter



Typeform



Facebook Pages



Google Drive



Asana



Dropbox

PREMIUM




Salesforce



RSS by Zapier



Webhooks by Zapier





Compare MANUAL vs. AUTOMATED qPCR-based detection and quantification of ZIKA virus using sciNote as ELN

MOTIVATION

Automate of qPCR-based ZIKV detection and quantification on two levels:


1. Using a robot (PIPETMAX® 268 automated pipetting workstation with qPCR assistant, Gilson Inc) & compared to manual pipetting
2. Using sciNote Open Source Electronic Lab Notebook (sciNote LLC) – to increase traceability




Trying to avoid as much paper as possible

Perform experiments in the lab that confirmed the link between ZIKV & Microcephaly

ZIKA VIRUS IS ASSOCIATED with MICROCEPHALY



The NEW ENGLAND JOURNAL of MEDICINE



Range of Microcephaly Severity



Baby with Typical Head Size

Baby with Microcephaly

Baby with Severe Microcephaly




BRIEF REPORT

Zika Virus Associated with Microcephaly

Jernej Mlakar, M.D., Misa Korva, Ph.D., Nataša Tul, M.D., Ph.D., Mara Popović, M.D., Ph.D., Mateja Poljšak-Prijatelj, Ph.D., Jerica Mraz, M.Sc., Marko Kolenc, M.Sc., Katarina Resman Rus, M.Sc., Tina Vesnaver Vipotnik, M.D., Vesna Fabjan Vodusek, M.D., Alenka Vizjak, Ph.D., Jože Pižem, M.D., Ph.D., Miroslav Petrovec, M.D., Ph.D., and Tatjana Avšič Županc, Ph.D.

EXPERIMENTAL DESIGN



SAMPLE LIST in sciNOTE



+ Add new sample Import Export

Show 25 entries Edit Delete Assign Unassign

Assigned	Sample name	Sample type	Sample group	Added on	Added by	Testing	Full description
<input checked="" type="checkbox"/>	Suriname	Cell culture	* RNA isolate	01.02.2017 17:39	Urška Čepin		Zika virus, strain Suriname
<input checked="" type="checkbox"/>	PRVABC59 COL	Cell culture	* RNA isolate	01.02.2017 17:39	Urška Čepin		Zika virus, strain PRVABC59 Puerto Rico
<input checked="" type="checkbox"/>	Pat neg	Patient	* RNA isolate	01.02.2017 17:39	Urška Čepin		Patient negative (blood)
<input checked="" type="checkbox"/>	Pat 5 semen	Patient	* RNA isolate	01.02.2017 17:39	Urška Čepin		Patient 5 (semen)
<input checked="" type="checkbox"/>	Pat 4 plasma	Patient	* RNA isolate	01.02.2017 17:39	Urška Čepin		Patient 2013 (plasma)
<input checked="" type="checkbox"/>	Pat 3 urine3	Patient	* RNA isolate	01.02.2017 17:39	Urška Čepin		Patient 3 (urine) - 3
<input checked="" type="checkbox"/>	Pat 3 urine2	Patient	* RNA isolate	01.02.2017 17:39	Urška Čepin		Patient 3 (urine) - 2
<input checked="" type="checkbox"/>	Pat 3 urine1	Patient	* RNA isolate	01.02.2017 17:39	Urška Čepin		Patient 3 (urine) - 1
<input checked="" type="checkbox"/>	Pat 2 urine	Patient	* RNA isolate	01.02.2017 17:39	Urška Čepin		Patient 2 (urine)
<input checked="" type="checkbox"/>	Pat 2 blood	Patient	* RNA isolate	01.02.2017 17:39	Urška Čepin		Patient 2 (blood)
<input checked="" type="checkbox"/>	Pat 1 brain	Patient	* RNA isolate	01.02.2017 17:39	Urška Čepin		Patient 1 (brain) 100x diluted
<input checked="" type="checkbox"/>	H/PF/2013 FP	Cell culture	* RNA isolate	01.02.2017 17:39	Urška Čepin		Zika virus, strain H/PF/2013 French Polynesia
<input checked="" type="checkbox"/>	FLR PRC	Cell culture	* RNA isolate	01.02.2017 17:39	Urška Čepin		Zika virus, strain FLR Columbia
<input checked="" type="checkbox"/>	976 Uganda	Cell culture	* RNA isolate	01.02.2017 17:39	Urška Čepin		Zika virus, strain 800/16 Brasil
<input checked="" type="checkbox"/>	800/16 Brasil	Cell culture	* RNA isolate	01.02.2017 17:39	Urška Čepin		Zika virus, strain 976 Uganda

EXAMPLE OF RNA ISOLATION PROTOCOL in sciNOTE

Start date: 15.02.2017 15:34 Due date: not set Status: In progress

Tags: Wet lab

EZ1 Virus Mini Kit v2.0 using EZ1 Advanced XL extraction system (Qiagen) is used for RNA from brain, plasma, blood, semen and cells culture samples.

(unlinked) [Load protocol](#) [Export protocol](#) [Save to protocol repository](#) [Complete task](#)

Protocol steps


[+ Add new](#) [Collapse all](#) [Expand all](#)

3 **Load cartridges, samples and labware into EZ1 instrument and START the instrument** | Published on 30.03.2017 17:06 by Matjaz Hren [↑](#) [↓](#) [↻](#) [🗑](#)

1. Open the instrument door.
2. Invert reagent cartridges 3 times to mix the magnetic particles. Then tap the cartridges to deposit the reagents to the bottom of their wells.
3. Follow the onscreen instructions for worktable setup, protocol variable selection, and data tracking:
 - **Note:** After sliding a reagent cartridge into the cartridge rack, press down on the cartridge until it clicks into place.
 - **Important:** make sure each elution tube for each sample tube is clearly labelled.
4. Close the instrument door.
5. Press "START" to start the protocol.

[Complete step](#)


Files




















qPCR ASSISTANT- AUTO GENERATED qPCR plate layout

	1	2	3	4	5	6	7	8	9	10	11	12
A	976 Uganda 3p	976 Uganda 3p	976 Uganda 3p	976 Uganda 3p	976 Uganda 3p	976 Uganda 3p	976 Uganda 3p	976 Uganda 3p	976 Uganda 3p	SC1	SC1	SC1
	1.0x	1.0x	1.0x	10.0x	10.0x	10.0x	100.0x	100.0x	100.0x	10 ⁵ x	10 ⁵ x	10 ⁵ x
B	800/16 Bra 3p	800/16 Bra 3p	800/16 Bra 3p	800/16 Bra 3p	800/16 Bra 3p	800/16 Bra 3p	800/16 Bra 3p	800/16 Bra 3p	800/16 Bra 3p	SC1	SC1	SC1
	1.0x	1.0x	1.0x	10.0x	10.0x	10.0x	100.0x	100.0x	100.0x	10 ⁵ x	10 ⁵ x	10 ⁵ x
C	FLR PRC 1p	FLR PRC 1p	FLR PRC 1p	FLR PRC 1p	FLR PRC 1p	FLR PRC 1p	FLR PRC 1p	FLR PRC 1p	FLR PRC 1p	Pat 2013 plaz	Pat 2013 plaz	Pat 3 uri3
	1.0x	1.0x	1.0x	10.0x	10.0x	10.0x	100.0x	100.0x	100.0x	1.0x	1.0x	1.0x
D	H/PP/2013 6p	H/PP/2013 6p	H/PP/2013 6p	H/PP/2013 6p	H/PP/2013 6p	H/PP/2013 6p	H/PP/2013 6p	H/PP/2013 6p	H/PP/2013 6p	Pat 1 brain	Pat 1 brain	Pat 5 spe
	1.0x	1.0x	1.0x	10.0x	10.0x	10.0x	100.0x	100.0x	100.0x	1.0x	1.0x	1.0x
E	PRVABC59 COL1p	PRVABC59 COL1p	PRVABC59 COL1p	PRVABC59 COL1p	PRVABC59 COL1p	PRVABC59 COL1p	PRVABC59 COL1p	PRVABC59 COL1p	PRVABC59 COL1p	Pat 2 blood	Pat 2 blood	Pat neg
	1.0x	1.0x	1.0x	10.0x	10.0x	10.0x	100.0x	100.0x	100.0x	1.0x	1.0x	1.0x
F	Suriname 5p	Suriname 5p	Suriname 5p	Suriname 5p	Suriname 5p	Suriname 5p	Suriname 5p	Suriname 5p	Suriname 5p	Pat 2 uri	Pat 2 uri	Pat 3 uri3
	1.0x	1.0x	1.0x	10.0x	10.0x	10.0x	100.0x	100.0x	100.0x	1.0x	1.0x	1.0x
G	NTC1	NTC1	NTC1	SC1	SC1	SC1	SC1	SC1	SC1	Pat 3 uri1	Pat 3 uri1	Pat 5 spe
	1.0x	1.0x	1.0x	10.0x	10.0x	10.0x	1000.0x	1000.0x	1000.0x	1.0x	1.0x	1.0x
H	SC1	SC1	SC1	SC1	SC1	SC1	SC1	SC1	SC1	Pat 3 uri2	Pat 3 uri2	Pat neg
	1.0x	1.0x	1.0x	100.0x	100.0x	100.0x	10 ⁴ x	10 ⁴ x	10 ⁴ x	1.0x	1.0x	1.0x

qPCR Assistant – layout of PIPETMAX tray for qPCR plate



qPCR Assistant – sample dilutions

	1	2	3	4	5	6	7	8	9	10	11	12
A	976 Uganda 3p 23,50 µL 1.0x 	Pat 3 uri3 17,60 µL 1.0x 	976 Uganda 3p 36,00 µL 10.0x	976 Uganda 3p 40,00 µL 100.0x	SC1ZIKV 1/10/100/ABI 7500 fast 1/10/100 7-log 36,00 µL 10 ⁵ x							
B	800/16 Bra 3p 23,50 µL 1.0x 	Pat 5 spe 17,60 µL 1.0x 	800/16 Bra 3p 36,00 µL 10.0x	800/16 Bra 3p 40,00 µL 100.0x	SC1ZIKV 1/10/100/ABI 7500 fast 1/10/100 7-log 40,00 µL 10 ⁶ x							
C	FLR PRC 1p 23,50 µL 1.0x 	Pat neg 17,60 µL 1.0x 	FLR PRC 1p 36,00 µL 10.0x	FLR PRC 1p 40,00 µL 100.0x	Pat 2013 plaz 17,60 µL 1.0x 							
D	H/PPF/2013 6p 23,50 µL 1.0x 		H/PPF/2013 6p 36,00 µL 10.0x	H/PPF/2013 6p 40,00 µL 100.0x	Pat 1 brain 17,60 µL 1.0x 							
E	PRVABC59 COL1p 23,50 µL 1.0x 		PRVABC59 COL1p 36,00 µL 10.0x	PRVABC59 COL1p 40,00 µL 100.0x	Pat 2 blood 17,60 µL 1.0x 							
F	Suriname 5p 23,50 µL 1.0x 		Suriname 5p 36,00 µL 10.0x	Suriname 5p 40,00 µL 100.0x	Pat 2 uri 17,60 µL 1.0x 							
G	NTC1ZIKV 1/10/100/ABI 7500 fast 1/10/100 7-log 23,10 µL 1.0x 		SC1ZIKV 1/10/100/ABI 7500 fast 1/10/100 7-log 36,00 µL 10.0x	SC1ZIKV 1/10/100/ABI 7500 fast 1/10/100 7-log 36,00 µL 1000.0x	Pat 3 uri1 17,60 µL 1.0x 							
H	SC1ZIKV 1/10/100/ABI 7500 fast 1/10/100 7-log 23,50 µL 1.0x 		SC1ZIKV 1/10/100/ABI 7500 fast 1/10/100 7-log 36,00 µL 100.0x	SC1ZIKV 1/10/100/ABI 7500 fast 1/10/100 7-log 36,00 µL 10 ⁴ x	Pat 3 uri2 17,60 µL 1.0x 							

qPCR Assistant – MASTER MIX RECIPES

Master Mix Recipes

Experiment template: 7500ABI-nondiluted_ABI7500Fast_96

Master mix 1: Assay: ZIKV 1

Required No. of reactions: 96

No. of reactions including extra volume: 110.2

No. of tubes:1

Reagent	Initial conc.	Final conc. / reaction	Number of reactions (µL)	
			1	110.2
water	-	-	9.2	1014.3
TaqMan Fast Virus 1-Step	4.0 x	1.0 x	5.0	551.2
BonnNS1 F	50.0 µM	0.63 µM	0.2	27.6
BonnNS1 R	50.0 µM	0.63 µM	0.2	27.6
BonnNS1 S	20.0 µM	0.3 µM	0.3	33.1
Total mix volume:			15.0 µL	1653.7 µL
Sample volume:			5.0 µL	
Total reaction volume:			20.0 µL	

qPCR Assistant – qPCR TEMPLATE

	1	2	3	4	5	6	7	8	9	10	11	12
A	976 Uganda 3p/1.0 U ZIKV 1/10/100	976 Uganda 3p/1.0 U ZIKV 1/10/100	976 Uganda 3p/1.0 U ZIKV 1/10/100	976 Uganda 3p/1.0 U ZIKV 1/10/100	976 Uganda 3p/1.0 U ZIKV 1/10/100	976 Uganda 3p/1.0 U ZIKV 1/10/100	976 Uganda 3p/1.0 U ZIKV 1/10/100	976 Uganda 3p/1.0 U ZIKV 1/10/100	976 Uganda 3p/1.0 U ZIKV 1/10/100	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0
B	800/16 Bra 3p/1.0 U ZIKV 1/10/100	800/16 Bra 3p/1.0 U ZIKV 1/10/100	800/16 Bra 3p/1.0 U ZIKV 1/10/100	800/16 Bra 3p/10.0 U ZIKV 1/10/100	800/16 Bra 3p/10.0 U ZIKV 1/10/100	800/16 Bra 3p/10.0 U ZIKV 1/10/100	800/16 Bra 3p/10.0 U ZIKV 1/10/100	800/16 Bra 3p/10.0 U ZIKV 1/10/100	800/16 Bra 3p/10.0 U ZIKV 1/10/100	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0
C	FLR PRC 1p/1.0 U ZIKV 1/10/100	FLR PRC 1p/1.0 U ZIKV 1/10/100	FLR PRC 1p/1.0 U ZIKV 1/10/100	FLR PRC 1p/10.0 U ZIKV 1/10/100	FLR PRC 1p/10.0 U ZIKV 1/10/100	FLR PRC 1p/10.0 U ZIKV 1/10/100	FLR PRC 1p/100.0 U ZIKV 1/10/100	FLR PRC 1p/100.0 U ZIKV 1/10/100	FLR PRC 1p/100.0 U ZIKV 1/10/100	Pat 2013 plaz/1.0 U ZIKV 1	Pat 2013 plaz/1.0 U ZIKV 1	Pat 3 uri3/1.0 U ZIKV 1
D	H/PF/2013 6p/1.0 U ZIKV 1/10/100	H/PF/2013 6p/1.0 U ZIKV 1/10/100	H/PF/2013 6p/1.0 U ZIKV 1/10/100	H/PF/2013 6p/10.0 U ZIKV 1/10/100	H/PF/2013 6p/10.0 U ZIKV 1/10/100	H/PF/2013 6p/10.0 U ZIKV 1/10/100	H/PF/2013 6p/10.0 U ZIKV 1/10/100	H/PF/2013 6p/10.0 U ZIKV 1/10/100	H/PF/2013 6p/10.0 U ZIKV 1/10/100	Pat 1 brain/1.0 U ZIKV 1	Pat 1 brain/1.0 U ZIKV 1	Pat 5 spe/1.0 U ZIKV 1
E	PRVABC59 COL1... U ZIKV 1/10/100	PRVABC59 COL1... U ZIKV 1/10/100	PRVABC59 COL1... U ZIKV 1/10/100	PRVABC59 COL1... U ZIKV 1/10/100	PRVABC59 COL1... U ZIKV 1/10/100	PRVABC59 COL1... U ZIKV 1/10/100	PRVABC59 COL1... U ZIKV 1/10/100	PRVABC59 COL1... U ZIKV 1/10/100	PRVABC59 COL1... U ZIKV 1/10/100	Pat 2 blood/1.0 U ZIKV 1	Pat 2 blood/1.0 U ZIKV 1	Pat neg/1.0 U ZIKV 1
F	Suriname 5p/1.0 U ZIKV 1/10/100	Suriname 5p/1.0 U ZIKV 1/10/100	Suriname 5p/1.0 U ZIKV 1/10/100	Suriname 5p/10.0 U ZIKV 1/10/100	Suriname 5p/10.0 U ZIKV 1/10/100	Suriname 5p/10.0 U ZIKV 1/10/100	Suriname 5p/100.0 U ZIKV 1/10/100	Suriname 5p/100.0 U ZIKV 1/10/100	Suriname 5p/100.0 U ZIKV 1/10/100	Pat 2 uri/1.0 U ZIKV 1	Pat 2 uri/1.0 U ZIKV 1	Pat 3 uri3/1.0 U ZIKV 1
G	NTC1 ZIKV 1/10/... N ZIKV 1/10/100	NTC1 ZIKV 1/10/... N ZIKV 1/10/100	NTC1 ZIKV 1/10/... N ZIKV 1/10/100	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0.1	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0.1	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0.1	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	Pat 3 uri1/1.0 U ZIKV 1	Pat 3 uri1/1.0 U ZIKV 1	Pat 5 spe/1.0 U ZIKV 1
H	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 1	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 1	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 1	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0.01	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0.01	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0.01	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	SC1 ZIKV 1/10/10... S ZIKV 1/10/100 0	Pat 3 uri2/1.0 U ZIKV 1	Pat 3 uri2/1.0 U ZIKV 1	Pat neg/1.0 U ZIKV 1

qPCR Assistant FILES UPLOADED TO sciNote

sciNOTE



Biosistemika



Hi, Matjaz Hren



Biosistemika / ZIKV application note / qPCR experiment: Robot... / Experiment setup with ...

PROTOCOLS

RESULTS

ACTIVITY

SAMPLES

REPORTS



Add new result: [Text](#) [Table](#) [File](#)

[Collapse all](#)

[Expand all](#)



qPCR plate pipetting sqLite protocol for PIPETMAX | Published on 30.03.2017 18:42 by Matjaz Hren



qPCR template for CYCLER | Published on 30.03.2017 18:41 by Matjaz Hren



qPCR pipetting guide | Published on 30.03.2017 18:40 by Matjaz Hren



Master mix sqLite protocol for PIPETMAX | Published on 30.03.2017 18:40 by Matjaz Hren



Master mix recipes | Published on 30.03.2017 18:39 by Matjaz Hren



Master mix pipetting guide | Published on 30.03.2017 18:39 by Matjaz Hren



Sample dilutions pipetting guide FOR MANUAL DILUTIONS | Published on 30.03.2017 18:38 by Matjaz Hren



Sample dilutions sqLite | Published on 30.03.2017 18:37 by Matjaz Hren




Sample dilutions pipetting guide | Published on 30.03.2017 18:36 by Matjaz Hren






RESULTS – MANUAL vs. ROBOT STANDARD CURVE


MANUAL



ROBOT



RESULTS – MANUAL vs. ROBOT ZIKV QUANTIFICATION




RESULTS – MANUAL vs. ROBOT ZIKV DETECTION

	Manual	PIPETMAX	Manual	PIPETMAX	Delta
Sample Name	Cq	Cq	AVG(Cq)	AVG(Cq)	Delta
Pat 1 brain/100.0	24.96	25.17	25.00	25.02	0.03
Pat 1 brain/100.0	25.03	24.88			
Pat 2 blood/1.0	32.62	34.45	32.78	34.00	1.22
Pat 2 blood/1.0	32.94	33.55			
Pat 2 urine/1.0	33.88	35.21	34.27	35.26	0.99
Pat 2 urine/1.0	34.66	35.31			
Pat 2013 plaz/1.0	33.55	34.58	33.55	36.15	2.60
Pat 2013 plaz/1.0	33.55	37.71			
Pat 3 urine1/1.0	34.06	36.12	34.98	36.94	1.97
Pat 3 urine1/1.0	35.90	37.77			
Pat 3 urine2/1.0	Na	Na		44.14	
Pat 3 urine2/1.0	Na	44.14			
Pat 3 urine3/1.0	Na	Na			
Pat 3 urine3/1.0	Na	Na			
Pat 5 semen/1.0	30.44	30.91	30.44	30.77	0.33
Pat 5 semen/1.0	30.44	30.63			
Pat negative/1.0	Na	Na			
Pat negative/1.0	Na	Na			




RESULTS – MANUAL vs. ROBOT SPEED

	Manual (min:s)	PIPETMAX (min:s)
Sample dilutions	21:30	11:30
Master mix	10:30	10:00
qPCR plate	32:00	16:30
SUM	64:00	38:00



SUMMARY

PIPETMAX qPCR Assistant



Patient samples

- Sample lists
- Patient data and sample tracking

Recipes and pipetting guide

- Auto-calculated volumes
- TRILUTION micro qPCR Assistant

qPCR reaction plate setup

- Automated with PIPETMAX® qPCR Assistant v1.3
- vs. manual liquid handling

Import run file into thermal cycler

- Run file generated by qPCR Assistant
- Reduces error vs. manual data entry

Real-time qPCR


- Applied Biosystems 7500
- Amplify with TaqMan® Gene Expression Assays

Data analysis

- Quality control checks
- Patient diagnosis



SUMMARY



SciNote Open source electronic lab notebook

CONCLUSIONS


using automation platform

- accurate standard dilution
- assay setup with significant increase in throughput

integration of data management software

- enables full traceability of samples and results → critical for accuracy in human diagnostics:
 - all in one place
 - all actions related to sample are recorded
 - backtrack any deviations/errors
- easier transfer of knowledge
- easier collaboration
- offers basic project management in science
- reduce possibility to lose data

OUTLOOK: BEYOND ROBOTS, BEYOND ELN



ACKNOWLEDGMENTS

Laboratory Automation and Data Management in Diagnostics – qPCR ZIKV Detection and Quantification

Klemen Zupančič¹, Miša Korva², Tatjana Avšič Županc², Urška Čepin¹,
Manca Pirc¹, Laura Simdon³, Matjaž Hren⁴

¹BioSistemika

²Institute of Microbiology and
Immunology, University of Ljubljana

³Gilson Inc.

⁴SciNote LLC

BioSistemika



Univerza v Ljubljani | Medicinska fakulteta

INŠTITUT ZA MIKROBIOLOGIJO IN IMUNOLOGIJO



GILSON[®]



SCINOTE

THANK YOU!



OPEN SOURCE ELECTRONIC LAB NOTEBOOK

www.scinote.net

GILSON/sciNote BOOTH #2

Matjaz Hren

E-mail: matjaz@scinote.net