



PIPETMAN® DIAMOND TIPS FOR PIPETMAX® 268 AND PIPETMAX® 278

AFFORDABLE AUTOMATION WITHOUT COMPROMISE

Unlike many automated liquid handlers that require proprietary consumables, Gilson's PIPETMAX® 268 and PIPETMAX® 278 automated liquid handlers use the same trusted **PIPETMAN® DIAMOND Tips** scientists rely on at the bench. You can automate critical workflows, like qPCR, normalization, or ELISA, without the added cost of special automation-only tips.

By supporting standard pipette tips, PIPETMAX® simplifies inventory management, reduces consumable costs, and makes automation more accessible. And if you already use other Gilson tools—such as PIPETMAN® pipettes—working within the same brand ecosystem helps reduce compatibility headaches, from purchasing and training to user interface familiarity.

ONE TIP FOR EVERY WORKFLOW

TIP TYPES:

- Standard, Filter, Sterile, 384-certified

PIPETTE HEADS:

- MAX8x20, MAX8x200, MAX4x1200

VOLUME RANGE:

- 1-1200 µL

CERTIFIED CLEAN TIPS:

- RNase/DNase-free
- Endotoxin-free sterile options

LOWER COST FOR LOWER IMPACT:

- One tip platform for manual and automated pipetting
- Reduce consumable costs and simplify inventory management



pipetman®

DIAMOND TIPS



pipetmax®

268



pipetmax®

278

DISPOSABLE TIP OPTIONS

QUANTITY: INCLUDES 10 PACKS OF 96 TIPS
MATERIAL: 100% VIRGIN POLYPROPYLENE

Part Number	Product Description	Tip Type	Volume Range
F167012	DL10 Reload Pack	Standard	0.1-20 µL
F167013	D200 Reload Pack	Standard	2-200 µL
F172200	DL10 Blister Refill	Standard	0.1-20 µL
F172300	D200 Blister Refill	Standard	2-200 µL
F172800	D1200 Blister Refill	Standard	50-1200 µL
F172801	D1200ST Blister Refill	Standard	50-1200 µL
F172203	DFL10ST Blister Refill	Sterile with filter	0.1-10 µL
F172303	DF30ST Blister Refill	Sterile with filter	2-30 µL
F172503	DF200ST Blister Refill	Sterile with filter	20-200 µL
F172803	DF1200ST Blister Refill	Sterile with filter	50-1200 µL
F172210	DSL10 Blister Refill	384 certified, sterile	0.1-20 µL
F172310	DS200 Blister Refill	384 certified, sterile	2-200 µL
F172211	DSL10ST Blister Refill	384 certified, sterile	0.1-20 µL
F172311	DS200ST Blister Refill	384 certified, sterile	2-200 µL
F172213	DSFL10ST Blister Refill	384 certified sterile with filter	0.1-10 µL
F172313	DSF30ST Blister Refill	384 certified sterile with filter	2-30 µL
F172513	DSF200ST Blister Refill	384 certified sterile with filter	20-200 µL

ORDERING INFORMATION

Visit [PIPETMAX Compatible Tips](#)